Protective Sealant for Treatment of Hypersensitive Dentin

TOKUYAMA SHIELD FORCE *PLUS*

Instructions for Use

Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <https://tokuyama-dental.com>

ENGLISH	2	БЪЛГАРСКИ	11	ČESKY	22	DEUTSCH	52
DANSK	43	EESTI KEEL	52	ESPAÑOL	62	SUOMI	72
FRANÇAIS	81	ΕΛΛΗΝΙΚΑ	91	HRVATSKI	103	MAGYAR	113
ITALIANO	123	LIETUVIŲ K.	133	LATVISKI	143	NEDERLANDS	152
NORSK	162	POLSKI	171	PORTUGUÊS	182	ROMÂNĂ	192
SVENSKA	202	SLOVENČINA	211	SLOVENŠČINA	221	TÜRKÇE	230

This product is certified as a medical device under the legislation of the European Union on medical devices by SGS CE1639, exclusively for the indication(s) shown in this IFU. Other non-medical uses ascribed to this product are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

Read all information, precautions and notes before using.

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

1. TOKUYAMA SHIELD FORCE PLUS is a one-component, single application, light-cured protective sealant designed for the treatment of hypersensitive dentin. TOKUYAMA SHIELD FORCE PLUS effectively infiltrates and seals exposed dentinal tubules to inhibit the flow of fluid through the tubules which is the cause of sensitivity. Coat of cured TOKUYAMA SHIELD FORCE PLUS is durable against mechanical wear, thermal stress, and chemical wear. TOKUYAMA SHIELD FORCE PLUS provides sufficient bond strength to teeth and resin-based materials such as composite resins, bonding agents, and multistep resin cements.
2. A light-curing unit with a wavelength range of 400 to 500 nm (peak: 470 nm) may be used for curing TOKUYAMA SHIELD FORCE PLUS.
3. TOKUYAMA SHIELD FORCE PLUS contains phosphoric acid monomer, bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), triethylene glycol dimethacrylate, 2-hydroxyethyl methacrylate (HEMA), camphorquinone, diphenyl(2,4,6-trimethylbenzoyl) phosphine oxide, alcohol and purified water.
4. TOKUYAMA SHIELD FORCE PLUS contains diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide as polymerization initiator in a concentration above 0.1% weight by weight. Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide is classified as reproductive toxicity 1B, in accordance with Part 3 of Annex VI to Regulation (EC) No 1272/2008 of the European Parliament and of the Council. Tokuyama Dental reviewed in detail the documents prepared by the regulatory authorities (Risk Assessment Committee of European Chemical Agency) and also conducted a risk assessment regarding TOKUYAMA SHIELD FORCE PLUS. It is concluded that TOKUYAMA SHIELD FORCE PLUS does not cause increased risk of adverse reproductive effects because the estimated exposures are below the tolerable daily intake (TDI values).
5. TOKUYAMA SHIELD FORCE PLUS is dispensed in a bottle.
6. Intended use environment: Dental clinics.

■ INTENDED PURPOSE

Intended purpose: Protective sealant for treatment of hypersensitive dentin

Patient target group: All age groups

■ CLINICAL BENEFIT

Relief of dentin hypersensitivity

■ INDICATIONS

TOKUYAMA SHIELD FORCE PLUS is indicated for:

- Treatment of hypersensitive dentin,
- Treatment of exposed cervical dentin caused through erosion or abrasion,
- Alleviation and/or prevention of tooth sensitivity after tooth preparation for direct and indirect restorations.

■ CONTRAINDICATIONS

- 1) TOKUYAMA SHIELD FORCE PLUS contains methacrylic monomers, organic solvents and acids. DO NOT use TOKUYAMA SHIELD FORCE PLUS for patients allergic to or hypersensitive to methacrylic monomers, related monomers, organic solvents and acids.
- 2) TOKUYAMA SHIELD FORCE PLUS is NOT indicated for:
 - Treatment of subgingival area where the gingival reattachment is expected,
 - Direct pulp capping,
 - Root canal sealing.
- 3) Not to be used if the patient suffers from severe gingival inflammation (especially conditions which can lead to tissue sloughing, e.g. acute necrotizing ulcerative gingivitis (ANUG)).

■ PRECAUTIONS

- 1) DO NOT use TOKUYAMA SHIELD FORCE PLUS for any purposes other than those listed herein. Use TOKUYAMA SHIELD FORCE PLUS only as directed.
- 2) TOKUYAMA SHIELD FORCE PLUS is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is suitable for use by non-dental care professionals.
- 3) DO NOT use TOKUYAMA SHIELD FORCE PLUS if the safety seals are broken or appear to have been tampered with.
- 4) If TOKUYAMA SHIELD FORCE PLUS causes an allergic reaction, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling TOKUYAMA SHIELD FORCE PLUS to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If TOKUYAMA SHIELD FORCE PLUS comes in contact with the examination gloves, remove

- and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 6) Avoid contact of TOKUYAMA SHIELD FORCE PLUS with eyes, mucous membrane, skin and clothing.
 - If TOKUYAMA SHIELD FORCE PLUS comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If TOKUYAMA SHIELD FORCE PLUS comes in contact with the mucous membrane, wipe the affected area immediately, and thoroughly flush the area with water. Affected areas may be whitened from protein coagulation, but such whitening tissue should disappear within 24 hours. If such whitening does not disappear after 24 hours, immediately contact a physician, and the patient should be so advised.
 - If TOKUYAMA SHIELD FORCE PLUS comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
 - Instruct the patient to rinse his/her mouth immediately after treatment.
 - 7) TOKUYAMA SHIELD FORCE PLUS should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
 - 8) To avoid the unintentional ingestion of TOKUYAMA SHIELD FORCE PLUS, do not leave it unsupervised within the reach of patients and children.
 - 9) DO NOT expose TOKUYAMA SHIELD FORCE PLUS or its vapor to open flame because TOKUYAMA SHIELD FORCE PLUS is flammable.
 - 10) To avoid cross infection, DO NOT reuse the disposable applicator included in the TOKUYAMA SHIELD FORCE PLUS package.
 - 11) Clean the dispensing well included in the TOKUYAMA SHIELD FORCE PLUS package thoroughly with alcohol after each use.
 - 12) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.

■ PRECAUTIONS FOR MEDICAMENTS AND MATERIALS

- 1) Some materials and medicaments (hemostatic solutions) inhibit curing of TOKUYAMA SHIELD FORCE PLUS for an extended period even after meticulous cleansing with water. DO NOT USE products which contain:
 - Eugenol
 - Hydrogen peroxide
 - Sodium hypochlorite
 - Diammine silver fluoride [molecular formula: $\text{Ag}(\text{NH}_3)_2\text{F}$]

- Phenols such as parachlorophenol, 2-methoxyphenol, phenol
 - Aluminum chloride
 - Ferric sulfate
 - Aluminum sulfate
 - Epinephrine
- 2) Application of phosphoric acid followed by application of TOKUYAMA SHIELD FORCE PLUS may reduce treatment efficacy.

■ STORAGE

- 1) Store under refrigerated conditions at temperatures between 0 to 10°C (32 to 50°F).
- 2) Keep away from heat, direct sunlight, sparks, and open flames.
- 3) DO NOT use after the indicated expiration date on the bottle/package.

■ DISPOSAL

Unused TOKUYAMA SHIELD FORCE PLUS is to be absorbed into an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

Refer to the packaging waste disposal guide shown on the end of this IFU.

■ CLINICAL PROCEDURES

1. Treatment of dentinal hypersensitivity

1.1. Cleaning

Clean tooth surface and remove plaque and saliva with a cotton pellet. A rubber cup with a fluoride-free paste under local anesthesia may also be used to clean.

1.2. Isolation

The rubber dam is the preferred method for isolation.

1.3. Drying

Dry the tooth by using a blotting technique or an air syringe technique.

- DO NOT desiccate the hypersensitive dentin. Desiccation can aggravate the hypersensitivity.
- Glass ionomer lining or calcium hydroxide should be applied if the area to be treated is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of TOKUYAMA SHIELD FORCE PLUS.

1.4. Dispensing

Open the bottle cap and dispense one or two drops of sealant into the Dispensing Well. Close the bottle cap tightly immediately after dispensing.

- Wipe off the excess sealant on the tip of nozzle before closing.

- Do not mix the sealant with other materials.

1.5. Application

Apply generous amounts of the sealant to the area to be treated using the disposable applicator. Leave undisturbed for 10 seconds or more. If the excess sealant pools on the area, blot the excess with a new disposable applicator.

- Protect the dispensed sealant and the inserted applicator from ambient light before the application using a light blocking plate.
- Complete the application within 5 minutes after dispensing because TOKUYAMA SHIELD FORCE PLUS contains a volatile alcohol.
- If saliva, blood, or other fluids contaminate the applied area, thoroughly rinse the area with water, dry and re-apply the sealant.
- Do not rinse the applied sealant with water except in unintentional contamination.

1.6. Air Dry

Using an oil-free air/water syringe, apply weak air flow continuously to the sealant surface until the runny sealant stays in the same position without any movement (usually for 5 seconds). Finish by using a strong air flow for 5 seconds or more. Use a vacuum aspirator to prevent spatter of the sealant.

- If spattering accidentally occurs, it may cause the tissue to whiten or a possible allergic reaction.
- To prevent spattering, refer to the following tips:
 - 1) To avoid inadvertent strong air;
 - a) Be sure to begin with a weak air flow outside the mouth,
 - b) Direct the air flow at the applied area.
 - 2) Extending the distance of the air/water syringe from the tooth reduces the air flow. The mirror reflection technique is also useful in reducing the air flow.

1.7. Light-cure

Keep the curing light tip within a distance of 2 mm from the surface. Then light-cure the surface for 10 seconds or more. If the surface is too large, divide the area into segments and light-cure each segment individually.

- Confirm that the light-curing unit has sufficient intensity ($>300 \text{ mW/cm}^2$) before using. Note that using a cracked light guide will lower the intensity.

1.8. Check for excess material

Check for excess material in the gingival sulcus, and remove the excess with a scaler or periodontal curette. Instruct the patient to rinse his/her mouth immediately after treatment.

1.9. Test for sensitivity

If a single application of TOKUYAMA SHIELD FORCE PLUS does not sufficiently relieve the hypersensitivity or sustained pain persists, the treatment can be repeated before alternative treatments are indicated. If hypersensitivity persists after a second application, seriously consider other sources of pain, such as traumatic occlusion, pulpal pathology.

- If additional restorations are required due to the esthetic reason, restore according to procedures in step “2.9. Adhesive treatment” and “2.10. Composite resin filling”.

1.10. Patient recall

Depending on individual oral hygiene, wear of the TOKUYAMA SHIELD FORCE PLUS may require re-application approximately every six months.

2. TOKUYAMA SHIELD FORCE PLUS sealant under direct restorations to alleviate and/or prevent tooth sensitivity after tooth preparation

2.1. Cleaning

Thoroughly clean the tooth surface using a rubber cup with a fluoride-free paste then rinse with water.

2.2. Isolation

The rubber dam is the preferred method for isolation.

2.3. Cavity preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in masking demarcations between the cavity margins and the restoration, enhancing both esthetics and retention.

2.4. Drying

Dry the cavity by using a blotting technique or an air syringe technique.

- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- Glass ionomer lining or calcium hydroxide should be applied if the area to be treated is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of TOKUYAMA SHIELD FORCE PLUS.

2.5. Dispensing

Open the bottle cap and dispense one or two drops of sealant into the Dispensing Well. Close the bottle cap tightly immediately after dispensing.

- Wipe off the excess sealant on the tip of nozzle before closing.
- Do not mix the sealant with other materials.

2.6. Application

Apply generous amounts of the sealant to the area to be treated using the disposable applicator. Leave undisturbed for 10 seconds or more. If the excess sealant pools on the area, blot the excess with a new disposable applicator.

- Protect the dispensed sealant and the inserted applicator from ambient light before the application using a light blocking plate.
- Complete the application within 5 minutes after dispensing because TOKUYAMA SHIELD FORCE PLUS contains a volatile alcohol.
- If saliva, blood, or other fluids contaminate the applied area, thoroughly rinse the area with water, dry and re-apply the sealant.
- Do not rinse the applied sealant with water except in unintentional contamination.

2.7. Air Dry

Using an oil-free air/water syringe, apply weak air flow continuously to the sealant surface until the runny sealant stays in the same position without any movement (usually for 5 seconds). Finish by using a strong air flow for 5 seconds or more. Use a vacuum aspirator to prevent spatter of the sealant.

- If spattering accidentally occurs, it may cause the tissue to whiten or a possible allergic reaction.
- To prevent spattering, refer to the following tips:
 - 1) To avoid inadvertent strong air;
 - a) Be sure to begin with a weak air flow outside the mouth,
 - b) Direct the air flow at the applied area.
 - 2) Extending the distance of the air/water syringe from the tooth reduces the air flow. The mirror reflection technique is also useful in reducing the air flow.

2.8. Light-cure

Keep the curing light tip within a distance of 2 mm from the surface. Then light-cure the surface for 10 seconds or more. If the surface is too large, divide the area into segments and light-cure each segment individually.

- Confirm that the light-curing unit has sufficient intensity ($>300 \text{ mW/cm}^2$) before using. Note that using a cracked light guide will lower the intensity.

2.9. Adhesive treatment

Immediately apply adhesive as described in manufacturer's directions for use.

- If using a total etch adhesive system, DO NOT etch the tooth surface prior to the use of TOKUYAMA SHIELD FORCE PLUS.

2.10. Composite resin filling

Place the restorative material over the cured adhesive according to manufacturer's directions for use.

3. TOKUYAMA SHIELD FORCE PLUS sealant under indirect restorations to alleviate and/or prevent tooth sensitivity after tooth preparation

3.1. Temporary Restoration

3.1.1. Tooth Preparation

Following tooth preparation, core build-up and fabrication of the provisional restoration, clean enamel and dentin with water spray and then air dry. Apply a rubber dam or other suitable isolation technique.

- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- Glass ionomer lining or calcium hydroxide should be applied if the area to be treated is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of TOKUYAMA SHIELD FORCE PLUS.

3.1.2. Dispensing

Open the bottle cap and dispense one or two drops of sealant into the Dispensing Well. Close the bottle cap tightly immediately after dispensing.

- Wipe off the excess sealant on the tip of nozzle before closing.
- Do not mix the sealant with other materials.

3.1.3. Application

Apply generous amounts of the sealant to the area to be treated using the disposable applicator. Leave undisturbed for 10 seconds or more. If the excess sealant pools on the area, blot the excess with a new disposable applicator.

- Protect the dispensed sealant and the inserted applicator from ambient light before the application using a light blocking plate.
- Complete the application within 5 minutes after dispensing because TOKUYAMA SHIELD FORCE PLUS contains a volatile alcohol.
- If saliva, blood, or other fluids contaminate the applied area, thoroughly rinse the area with water, dry and re-apply the sealant.
- Do not rinse the applied sealant with water except in unintentional contamination.

3.1.4. Air Dry

Using an oil-free air/water syringe, apply weak air flow continuously to the sealant surface until the runny sealant stays in the same position without any movement (usually for 5 seconds). Finish by using a strong air flow for 5 seconds or more. Use a vacuum aspirator to

prevent spatter of the sealant.

- If spattering accidentally occurs, it may cause the tissue to whiten or a possible allergic reaction.
- To prevent spattering, refer to the following tips:
 - 1) To avoid inadvertent strong air;
 - a) Be sure to begin with a weak air flow outside the mouth,
 - b) Direct the air flow at the applied area.
 - 2) Extending the distance of the air/water syringe from the tooth reduces the air flow. The mirror reflection technique is also useful in reducing the air flow.

3.1.5. Light-cure

Keep the curing light tip within a distance of 2 mm from the surface. Then light-cure the surface for 10 seconds or more. If the surface is too large, divide the area into segments and light-cure each segment individually.

- Confirm that the light-curing unit has sufficient intensity ($>300 \text{ mW/cm}^2$) before using. Note that using a cracked light guide will lower the intensity.

3.1.6. Dismissal

Wipe the oxygen inhibition layer with alcohol-soaked cotton swab or gauze.

3.1.7. Impression Taking

Take impression with impression material according to manufacturer's instructions.

3.1.8. Temporary Cementing

Seat the temporary/provisional restoration on the tooth with temporary cement (water-based cement) according to manufacturer's instructions. Temporary stopping (gutta-percha) may also be used for temporary filling.

- Do not use resin-based material (cement/filling) because it may adhere to the cured TOKUYAMA SHIELD FORCE PLUS.

3.2. Permanent Restoration

3.2.1. Cleaning of Cavity/Abutment

Remove temporary restoration, residual temporary cement and hard deposits from the tooth by ultrasonic scaler or other cleaning instruments. Thoroughly clean all contaminants from the tooth surface to be bonded. Rinse with water and dry in the usual manner.

- Elements listed below, which inhibit the setting of resin cement, should be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol, citric acid, EDTA, or the application of phosphoric acid for 2 to 3 seconds before cementation.

- 1) Silicone oil from fit check material

- 2) Oil mist from handpiece
- 3) Saliva, blood and exudates

3.2.2. Cementing

Cement the restoration with resin cement according to its manufacturer's instructions.

- The use of multistep resin cements are recommended for cementing. The use of self-adhesive resin cements are not recommended due to potential lowered adhesion of the cement.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The manufacturer of the TOKUYAMA SHIELD FORCE PLUS is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.

Product specifications of the TOKUYAMA SHIELD FORCE PLUS are subject to change without notification. When the product specifications are changed, the instructions and precautions may also be changed.

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Packaging waste disposal guide

(According to European Commission Decision 97/129/EC)

ENGLISH

Packaging waste disposal guide

Bottle (cap and body)	Dispose of as industrial waste.
Case body for applicators	CP7 (Cellulose propionate): Recycling. Follow local instructions for disposal.
Case cap for applicators	PP5 (Polypropylene): Recycling. Follow local instructions for disposal.
Paperbox and/or partation	PAP21 (Non-corrugated fibreboard): Recycling. Follow local instructions for disposal.
IFU and/or guide	PAP22 (Paper): Recycling. Follow local instructions for disposal.

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