Resin-based Dental Restorative Material

ESTELITE BULK FILL *Flow* SYRINGE**ESTELITE BULK FILL *Flow*** PLT

Instructions for Use

**Tokuyama Dental Corporation**

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <https://tokuyama-dental.com>

ENGLISH	2	БЪЛГАРСКИ	5	ČESKY	9	DEUTSCH	13
DANSK	17	EESTI KEEL	21	ESPAÑOL	25	SUOMI	29
FRANÇAIS	32	ΕΛΛΗΝΙΚΑ	37	HRVATSKI	41	MAGYAR	45
ITALIANO	49	LIETUVIŲ K.	53	LATVISKI	56	NEDERLANDS	60
NORSK	64	POLSKI	68	PORTUGUÊS	72	ROMÂNĂ	76
SVENSKA	80	SLOVENČINA	83	SLOVENŠČINA	87	TÜRKÇE	91

This product is certified as a medical device under the legislation of the European Union on medical devices by SGS CE1639, exclusively for the indication(s) shown in this IFU. Other non-medical uses ascribed to this product are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

Read all information, precautions and notes before using.

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

- 1) ESTELITE BULK FILL Flow is a low viscosity, light-cured, radiopaque composite resin. This low stress flowable material can be placed 4mm increments. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049.
- 2) ESTELITE BULK FILL Flow incorporates Radical-Amplified Photopolymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time.
- 3) ESTELITE BULK FILL Flow contains 70% by weight (56% by volume) of spherical silica-zirconia filler and composite filler. The composite filler offers reduced polymerization shrinkage stress due to its round shape. The spherical filler contained in ESTELITE BULK FILL Flow (mean particle size: 200nm, particle size range: 100 to 300nm) facilitates excellent gloss retention and wear resistance. ESTELITE BULK FILL Flow contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate (TEGDMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), Mequinol, Dibutyl hydroxyl toluene and UV absorber.
- 4) ESTELITE BULK FILL Flow is provided either in SYRINGE or Pre Loaded Tip (PLT).

■ SHADE

U (Universal), A1, A2, A3, B1

- Please check with your local distributor for available shades.

■ INDICATIONS

- Direct anterior and posterior restorations
- Cavity lining
- Blocking out cavity undercuts before fabrication indirect restorations
- Repair of porcelain/composite

■ CONTRAINDICATIONS

ESTELITE BULK FILL Flow contains methacrylic monomers. DO NOT use ESTELITE BULK FILL Flow for patients allergic to or hypersensitive to methacrylic and related monomers or to any of the other ingredients.

■ PRECAUTIONS

- 1) DO NOT use ESTELITE BULK FILL Flow for any purposes other than those listed in these instructions. Use ESTELITE BULK FILL Flow only as directed herein.
- 2) ESTELITE BULK FILL Flow is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
- 3) DO NOT use ESTELITE BULK FILL Flow if the safety seals are broken or appear to have been tampered with.
- 4) If ESTELITE BULK FILL Flow causes an allergic reaction or hypersensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl, or latex) at all times when handling ESTELITE BULK FILL Flow to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE BULK FILL Flow comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 6) Avoid contact of ESTELITE BULK FILL Flow with eyes, mucosal membrane, skin and

clothing.

- If ESTELITE BULK FILL Flow comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE BULK FILL Flow comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If ESTELITE BULK FILL Flow comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
- 7) ESTELITE BULK FILL Flow should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
 - 8) To avoid the unintentional ingestion of ESTELITE BULK FILL Flow, do not leave it unsupervised within the reach of patients and children.
 - 9) Clean the placement instruments and brushes with alcohol after use.
 - 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
 - 11) ESTELITE BULK FILL Flow is available in semi-translucent shades. If applied in cavity Class III or IV direct restorations, the cured ESTELITE BULK FILL Flow may appear dark. Avoid using ESTELITE BULK FILL Flow for the direct restoration of Class III or IV when opacity is paramount.
 - 12) Light-cure ESTELITE BULK FILL Flow extra orally to confirm the time needed for complete hardening of ESTELITE BULK FILL Flow with your light-curing unit prior to doing the restoration.
 - 13) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

■ CURING TIME AND DEPTH

ESTELITE BULK FILL Flow enables a 4mm depth of cure (U shade) by light-curing for 10 seconds with Halogen or LED light-curing unit (Light intensity: 800mW/cm² or more, wavelength: 400-500nm). The depth of cure was determined based on test results performed in accordance with section 7.10 “depth of cure” of ISO 4049.

■ SPECIAL NOTES FOR THE USE OF PLT

- 1) PLTs are designed for single patient use only. To avoid cross infection, Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.
- 2) Dispensers are not provided for ESTELITE BULK FILL Flow. Use a dispenser that fits the PLT of ESTELITE BULK FILL Flow. For proper use and proper disinfection, see manufacturer’s instructions.
- 3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

■ STORAGE

- 1) Store ESTELITE BULK FILL Flow at a temperature between 0 - 25°C (32 - 77°F).
- 2) AVOID direct exposure to light and heat.
- 3) DO NOT use ESTELITE BULK FILL Flow after the indicated date of expiration on the syringe or PLT package.

■ DISPOSAL

To safely dispose of excess ESTELITE BULK FILL Flow, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

Do not allow to enter drains, sewers or watercourses.

Must not be disposed of together with household garbage.

Follow local instructions for disposal.

Refer to the packaging waste disposal guide shown on the end of this IFU.

■ CLINICAL PROCEDURE

1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Shade Selection

Prior to isolation of tooth, select the appropriate shade of ESTELITE BULK FILL Flow using shade card provided for ESTELITE BULK FILL Flow.

3. Isolation

A rubber dam is the preferred method of isolation.

4. Cavity Preparation

Prepare the cavity and rinse with water.

- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

5. Pulp Protection

Calcium hydroxide should be applied if the cavity is in close proximity to the pulp. **DO NOT USE EUGENOL-BASED MATERIALS** to protect the pulp as these materials will inhibit curing ESTELITE BULK FILL Flow. ESTELITE BULK FILL Flow is not indicated for direct pulp capping.

6. Bonding system

Apply a bonding system according to its manufacturer's instructions.

7. Dispensing


7-1. PLT

Please read ■ SPECIAL NOTES FOR THE USE OF PLT before use.

- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.

7-2. SYRINGE

- 1) Hold the nozzle and turn the cap counterclockwise for removal. [1]
- 2) Hold the nozzle and place the Dispensing Tip onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

- Accessories and Dispensing Tip are exclusively designed for ESTELITE BULK FILL Flow. Avoid using Dispensing Tips from other brands with ESTELITE BULK FILL Flow.
- Dispensing Tips are also sold separately.
- Use gloves when attaching the Dispensing Tip.

8. Filling and contouring

Cavity can be directly filled from the Dispensing Tip/PLT. Increments should not exceed 4mm at a time.

- **DO NOT** apply excessive force to extrude the material from the syringe or PLT immediately after removal from the refrigerator.

- DO NOT mix ESTELITE BULK FILL Flow with other brands of resin composite or other shades of ESTELITE BULK FILL Flow to avoid incomplete cure or entrapment of air bubbles.
- In case of syringe, to avoid cross infection, after extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.

9. Curing

Light-cure each increment for at least the indicated time (please see the curing recommendations below). Light-cure each layer (up to 4mm increment) using Halogen light-curing or LED light-curing unit with a wavelength of 400-500 nm. Be sure the light-curing unit with an intensity of 800mW/cm² or more. Hold the light guide tip as close to the restorative as possible during light exposure.

CURING RECOMMENDATIONS

Shade	Increment Depth	Light Output	Curing Time
All shades	4mm	≥800mW/cm ²	10 sec.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

10. Finishing

Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

- Instruct the patient adequately on the maintenance of the restoration after treatment.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.

Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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Packaging waste disposal guide

(According to European Commission Decision 97/129/EC)

ENGLISH

Packaging waste disposal guide

PLT (cap and body)	Dispose of as industrial waste.
SYRINGE (cap and body)	Dispose of as industrial waste.
Packaging bag for PLT	PP5 (Polypropylene): Recycling. Follow local instructions for disposal.
Packaging bag for tip	LDPE4 (Low density polyethylene): Recycling. Follow local instructions for disposal.
Paperbox and/or partation	PAP21 (Non-corrugated fibreboard): Recycling. Follow local instructions for disposal.
IFU and/or guide	PAP22 (Paper): Recycling. Follow local instructions for disposal.

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