Light-cured Characterizing Material

ESTELITE COLOR

CE 1639

Instructions for Use

Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <https://tokuyama-dental.com>

ENGLISH	2	БЪЛГАРСКИ	5	ČESKY	9	DEUTSCH	13
DANSK	17	EESTI KEEL	20	ESPAÑOL	24	SUOMI	28
FRANÇAIS	31	ΕΛΛΗΝΙΚΑ	35	HRVATSKI	39	MAGYAR	43
ITALIANO	47	LIETUVIŲ K.	51	LATVISKI	54	NEDERLANDS	58
NORSK	62	POLSKI	65	PORTUGUÊS	69	ROMÂNĂ	73
SVENSKA	77	SLOVENČINA	80	SLOVENŠČINA	84	TÜRKÇE	88

This product is certified as a medical device under the legislation of the European Union on medical devices by SGS CE1639, exclusively for the indication(s) shown in this IFU. Other non-medical uses ascribed to this product are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

Carefully read all instructions, precautions and notes before using ESTELITE COLOR.

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

- 1) ESTELITE COLOR is a low viscosity, light-cured tint for individualized shade characterizations of direct and indirect resin restorations. ESTELITE COLOR is applied below or in-between the layers of the composite restoration.
- 2) ESTELITE COLOR contains 58% by weight of silica-zirconia filler and silica-titania filler. The monomer matrix contains Bisphenol A polyethoxy methacrylate (Bis-MPEPP), Triethylene glycol dimethacrylate (TEGDMA), and 1,6-bis(methacryl-ethyloxycarbonylamino) trimethyl hexane (UDMA).

■ SHADE

ESTELITE COLOR contains 9 Tints and 4 Opaquers.

Tints: Blue, Gray, Lavender, Ochre, Dark Brown, White, Yellow, Red, Clear

Opaquers: Pink Opaque, Low Chroma Opaque, Medium Chroma Opaque, High Chroma Opaque

Low Chroma Opaque (LCO): for shade range A1, B1, C1 and Bleached teeth

Medium Chroma Opaque (MCO): for shade range A2, B2, D2, and C2

High Chroma Opaque (HCO): for shade range A3, B3, D3, C3 and darker

■ INDICATIONS

- Individualized shade characterizations of direct and indirect resin restorations.
- Masking of discoloration.

■ CONTRAINDICATIONS

ESTELITE COLOR contains methacrylic monomers. DO NOT use ESTELITE COLOR for patients allergic to or hypersensitive to methacrylic and related monomers.

■ PRECAUTIONS

- 1) DO NOT use ESTELITE COLOR for any purposes other than those listed in these instructions. Use ESTELITE COLOR only as directed herein.
- 2) ESTELITE COLOR is designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
- 3) DO NOT use ESTELITE COLOR if the safety seals are broken or appear to have been tampered with.
- 4) If ESTELITE COLOR causes an allergic reaction or hypersensitivity, discontinue use immediately.
- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling ESTELITE COLOR to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE COLOR comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 6) Avoid contact of ESTELITE COLOR with eyes, mucosal membrane, skin and clothing.
 - If ESTELITE COLOR comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE COLOR comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If ESTELITE COLOR comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
- 7) ESTELITE COLOR should not be ingested or aspirated. Ingestion or aspiration may cause

serious injury.

- 8) To avoid the unintentional ingestion of ESTELITE COLOR, do not leave it unsupervised within the reach of patients and children.
- 9) Clean the placement instruments with alcohol after use.
- 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.

■ STORAGE

- 1) Store ESTELITE COLOR at a temperature between 0 - 25°C (32 - 77°F).
- 2) AVOID direct exposure to light and heat.
- 3) DO NOT use ESTELITE COLOR after the indicated date of expiration on the syringe/package.

■ DISPOSAL

To safely dispose of excess ESTELITE COLOR, extrude unused portion from the syringe and light-cure before disposal.

Do not allow to enter drains, sewers or watercourses.


Must not be disposed of together with household garbage.

Follow local instructions for disposal.

Refer to the packaging waste disposal guide shown on the end of this IFU.

■ ATTACHING THE SYRINGE TIP

- 1) Hold the nozzle and turn the cap counterclockwise for removal. [1]
 - 2) Hold the nozzle and place the Syringe Tip onto it. Turn the Syringe Tip clockwise until it is securely locked. [2]
- Accessories and Syringe Tip are designed for ESTELITE COLOR and PALFIQUE ESTELITE LV. Avoid using Syringe Tips from other brands with ESTELITE COLOR.
 - Syringe Tips are also sold separately.
 - Use gloves when attaching the Syringe Tip.

■ CLINICAL PROCEDURE

Individualized shade characterizations of direct resin restorations

1. Cavity preparation & Bonding system

Prepare the cavity and rinse with water. Apply light-cured bonding system according to manufacturer's instructions.

2. Dentin shade resin composite

Based on the multilayer approach, a dentin shade resin composite is applied over the bonded surface. If required, the mamelons or fissures are sculpted with a fine tipped instrument at this stage. Then the dentin is light cured according to manufacturer's instructions.

3. Dispensing the ESTELITE COLOR

Extrude the small amount of ESTELITE COLOR on the Disposable Dish. If required, mix ESTELITE COLOR with other shades of ESTELITE COLOR. After extruding the material, wipe the nozzle with gauze and replace the cap immediately.

- Tips for shade selection
 - Incisal edge: Blue (youthful look), Gray (mature look), Lavender (elderly look)

- Increasing saturation: Ochre (fossae areas), Dark Brown (pits and fissures; stained craze lines)
- White spots: White
- Color adjustment: Yellow, Ochre, Red (mixed in varying proportions)
- Chroma adjustment: Clear (mixed with another shade)
- Neutralizing grayness: Pink Opaque (mixed with LCO, MCO, or HCO)

Tints can be mixed in varying proportions to achieve the correct hue, value and chroma, as needed. Opaquers can be mixed in varying proportions to achieve the correct value and chroma, as needed.

- Protect the dispensed material from ambient light before the application using a Light Shield Cover.
- DO NOT mix ESTELITE COLOR with other brands of products to avoid incomplete cure.

4. Applying the ESTELITE COLOR

Using a fine tipped instrument or a fine tipped brush, apply ESTELITE COLOR sparingly. The maximum thickness is 0.5 mm (Dark Brown: 0.05 mm).

5. Curing

Use a light curing unit with a wavelength range of 400 to 500 nm (peak: 470 nm, intensity $\geq 500 \text{ mW/cm}^2$). Light-cure for at least 30 seconds, ensure the curing light tip does not exceed 2 mm from the surface.

6. Enamel shade resin composite

An enamel shade resin composite is applied over the colored surface. Then the enamel is light cured; finish and polish according to manufacturer's instructions. ESTELITE COLOR must be placed below or in-between the layers of the composite restoration.

- Instruct the patient adequately on the maintenance of the restoration after treatment.

Masking of discoloration

1. Cavity preparation & Bonding system

Prepare the cavity and rinse with water. Apply light-cured bonding system according to manufacturer's instructions.

2. Dispensing the ESTELITE COLOR

Extrude the small amount of ESTELITE COLOR on the Disposable Dish. If required, mix ESTELITE COLOR with other shades of ESTELITE COLOR. After extruding the material, wipe the nozzle with gauze and replace the cap immediately.

- Tips for shade selection
 - Masking dark areas: LCO, MCO, HCO
- Protect the dispensed material from ambient light before the application using a Light Shield Cover.
- DO NOT mix ESTELITE COLOR with other brands of products to avoid incomplete cure.
- DO NOT use self-cured bonding systems. We advise that light-cured bonding systems should be used with ESTELITE COLOR since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

3. Applying ESTELITE COLOR

Using a fine tipped brush or other suitable instruments, apply the appropriate Opaquer of ESTELITE COLOR evenly over the discolored areas to elevate the color value. The maximum thickness is 0.5 mm for each increment.

- DO NOT apply ESTELITE COLOR directly onto the tooth without application of bonding agent.
- Tips for thickness of Opaquer

The required average thickness of Opaquer is 0.35 mm to mask the discolored areas

adequately.

4. Curing

Use a light curing unit with a wavelength range of 400 to 500 nm (peak: 470 nm, intensity $\geq 500 \text{ mW/cm}^2$). Light-cure for at least 30 seconds, ensure the curing light tip does not exceed 2 mm from the surface.

5. Resin composite

A resin composite is applied over the material. Then the composite is light cured; finish and polish according to manufacturer's instructions. ESTELITE COLOR must be placed below or in between the layers of the composite restoration.

Individualized shade characterizations of indirect resin restorations

Use indirect composite resin system according to manufacture's instructions. During the fabrication, ESTELITE COLOR is placed between the layers of resin-based restorative material according to the procedure described for direct restorations (direct staining technique).

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.

Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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Packaging waste disposal guide

(According to European Commission Decision 97/129/EC)

ENGLISH

Packaging waste disposal guide

SYRINGE (cap and body)	Dispose of as industrial waste.
Packaging bag for tip, cap and Disposable Dish	LDPE4 (Low density polyethylene): Recycling. Follow local instructions for disposal.
Paperbox and partation	PAP21 (Non-corrugated fibreboard): Recycling. Follow local instructions for disposal.
IFU and guide	PAP22 (Paper): Recycling. Follow local instructions for disposal.

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