Resin-based Dental Restorative Material

**OMNICHROMA *Flow* BULK**

SYRINGE

**Instructions for Use****Tokuyama Dental Corporation**

38-9, Taitou 1-chome, Taitou-ku, Tokyo, Japan

Tel: +81-3-3835-7201

URL: <https://tokuyama-dental.com>

|                 |           |                    |           |                    |           |                   |            |
|-----------------|-----------|--------------------|-----------|--------------------|-----------|-------------------|------------|
| <b>ENGLISH</b>  | <b>2</b>  | <b>БЪЛГАРСКИ</b>   | <b>6</b>  | <b>ČESKY</b>       | <b>11</b> | <b>DEUTSCH</b>    | <b>15</b>  |
| <b>DANSK</b>    | <b>19</b> | <b>EESTI KEEL</b>  | <b>24</b> | <b>ESPAÑOL</b>     | <b>28</b> | <b>SUOMI</b>      | <b>32</b>  |
| <b>FRANÇAIS</b> | <b>36</b> | <b>ΕΛΛΗΝΙΚΑ</b>    | <b>41</b> | <b>HRVATSKI</b>    | <b>46</b> | <b>MAGYAR</b>     | <b>50</b>  |
| <b>ITALIANO</b> | <b>55</b> | <b>LIETUVIŲ K.</b> | <b>59</b> | <b>LATVISKI</b>    | <b>63</b> | <b>NEDERLANDS</b> | <b>67</b>  |
| <b>NORSK</b>    | <b>72</b> | <b>POLSKI</b>      | <b>76</b> | <b>PORTUGUÊS</b>   | <b>81</b> | <b>ROMÂNĂ</b>     | <b>85</b>  |
| <b>SVENSKA</b>  | <b>90</b> | <b>SLOVENČINA</b>  | <b>94</b> | <b>SLOVENŠČINA</b> | <b>98</b> | <b>TÜRKÇE</b>     | <b>102</b> |

This product is certified as a medical device under the legislation of the European Union on medical devices by SGS CE1639, exclusively for the indication(s) shown in this IFU. Other non-medical uses ascribed to this product are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

Carefully read all information, precautions, notes and materials before using.

## ■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

- 1) OMNICHROMA FLOW BULK is a low viscosity, light-cured, radiopaque composite resin for use in anterior and posterior restorations and is indicated for all carious classes. OMNICHROMA FLOW BULK is a single shade material. Direct placement of OMNICHROMA FLOW BULK into a prepared cavity allows for easy handling. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049.
- 2) OMNICHROMA FLOW BULK is a flowable material which can be placed 3.5mm increments.
- 3) OMNICHROMA FLOW BULK contains 69% by weight (55% by volume) of spherical silica-zirconia filler (mean particle size: 0.3µm, particle size range: 0.2 to 0.4µm) and composite filler. OMNICHROMA FLOW BULK contains 1,6-bis(methacryl-ethyloxycarbonylamino) trimethyl hexane (UDMA), Triethylene glycol dimethacrylate (TEGDMA), Mequinal, Dibutyl hydroxyl toluene and UV absorber.
- 4) OMNICHROMA FLOW BULK is available in SYRINGE.

## ■ INDICATIONS

- Direct anterior and posterior restorations
- Cavity base or liner
- Blocking out cavity undercuts before fabricating indirect restorations
- Repair of porcelain/composite

## ■ CONTRAINDICATIONS

OMNICHROMA FLOW BULK contains methacrylic monomers and UV absorber. DO NOT use OMNICHROMA FLOW BULK for patients allergic to or hypersensitive to methacrylic and related monomers, UV absorber or any of the other ingredients.

## ■ PRECAUTIONS

- 1) DO NOT use OMNICHROMA FLOW BULK for any purpose other than those listed in these instructions. Use OMNICHROMA FLOW BULK only as directed herein.
- 2) OMNICHROMA FLOW BULK is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
- 3) DO NOT use OMNICHROMA FLOW BULK if the safety seals are broken or appear to have been tampered with.
- 4) If OMNICHROMA FLOW BULK causes an allergic reaction or hypersensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling OMNICHROMA FLOW BULK to avoid the possibility of allergic reaction from methacrylic monomers or UV absorber.

Note: Certain substance/materials may penetrate through examination gloves. If OMNICHROMA FLOW BULK comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.

- 6) Avoid contact of OMNICHROMA FLOW BULK with eyes, mucosal membrane, skin and clothing.
  - If OMNICHROMA FLOW BULK comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
  - If OMNICHROMA FLOW BULK comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.

- If OMNICHROMA FLOW BULK comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
  - Instruct the patient to rinse his mouth immediately after treatment.
- 7) OMNICHROMA FLOW BULK should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
  - 8) To avoid the unintentional ingestion of OMNICHROMA FLOW BULK, do not leave it unsupervised within the reach of patients and children.
  - 9) Clean the placement instruments and brushes with alcohol after use.
  - 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
  - 11) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude. After removing from the refrigerator, allow the product to sit 15 minutes or until it reaches room temperature.
  - 12) Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.

## ■ INDICATIONS FOR FILLING AND CURING

OMNICHROMA FLOW BULK is designed to be cured by either a halogen or LED curing-light with a wavelength of 400-500 nm. OMNICHROMA FLOW BULK enables a 3.5mm increment depth by light-curing for 20 seconds with Halogen or LED light-curing unit (Light intensity: 800mW/cm<sup>2</sup> or more, wavelength: 400-500nm).

Be sure to light-cure OMNICHROMA FLOW BULK extra-orally and check the time needed for complete hardening of OMNICHROMA FLOW BULK with your light-curing unit before performing the bonding procedure. The following table summarizes the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

| Light type | Intensity (mW/cm <sup>2</sup> ) | Curing time (seconds) | Increment depth (mm) <sup>1)</sup> |
|------------|---------------------------------|-----------------------|------------------------------------|
| Halogen    | 600                             | 30                    | 3.7                                |
|            | 800                             | 20                    | 3.6                                |
|            |                                 | 30                    | 4.1                                |
| LED        | 800                             | 20                    | 3.5                                |
|            |                                 | 30                    | 3.9                                |
|            | 900                             | 20                    | 3.5                                |
|            |                                 | 30                    | 4.0                                |
|            | 1000                            | 20                    | 3.6                                |
|            |                                 | 30                    | 4.1                                |
|            | 1200                            | 20                    | 3.7                                |
|            |                                 | 30                    | 4.2                                |

- 1) Increment depth was determined on the basis of test results performed in accordance with “depth of cure” of ISO 4049.

## ■ STORAGE

- 1) Store OMNICHROMA FLOW BULK at temperatures between 0 - 25°C (32 - 77°F).
- 2) AVOID direct exposure to light and heat.
- 3) DO NOT use OMNICHROMA FLOW BULK after the indicated date of expiration on the syringe package.

## ■ DISPOSAL

To safely dispose of excess OMNICHROMA FLOW BULK, extrude unused portion from SYRINGE and light-cure before disposal.

Do not allow to enter drains, sewers or watercourses.

Must not be disposed of together with household garbage.

Follow local instructions for disposal.

Refer to the packaging waste disposal guide shown on the end of this IFU.

## ■ CLINICAL PROCEDURE

Determine if this product is suitable for application to the case before use. The desirable esthetic outcome may not be obtained depending on the case (e.g. a case where very little tooth structure remains, the underlying tooth is extremely discolored, or the color of the underlying tooth is to be changed to another color, etc.).

### 1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

### 2. Isolation

A rubber dam is the preferred method of isolation.

### 3. Cavity Preparation

Prepare the cavity. In case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.

Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II), bevels and chamfers assist in erasing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.

- Scalloped bevel could be desirable depending on the esthetics and retention.

- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent, followed by a metal primer if metal is exposed on the fractured surfaces according to manufacturer's instructions.

### 4. Pulp Protection

Calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing OMNICHROMA FLOW BULK. OMNICHROMA FLOW BULK is not indicated for direct pulp capping.

### 5. Bonding System

Apply a bonding system according to its manufacturer's instructions.


- If using self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

- In case of porcelain/composite repairs, be sure to condition the surface of porcelain/composite with a silane coupling agent BEFORE performing the bonding procedure.

### 6. Dispensing

- Hold the nozzle and turn the cap counterclockwise for removal. [1]

- Hold the nozzle and place the Dispensing Tip (metal tip) onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

- Accessories and Dispensing Tips are exclusively designed for the SYRINGE of OMNICHROMA FLOW BULK. Use the attached Dispensing Tips (metal tips).
- Dispensing Tips (TOKUYAMA Dispensing Tip) are also sold separately.
- Use uncontaminated gloves when handling the syringe and attaching the Dispensing Tip.

## 7. Filling and contouring

Cavities can be directly filled from the Dispensing Tip, or indirectly filled with an instrument after extruding the paste on the mixing pad. Increments should not exceed the indicated curing depth at a time (please refer to the aforementioned table).

- DO NOT apply excessive force to extrude the material from the syringe immediately after removal from the refrigerator.
- DO NOT mix OMNICHROMA FLOW BULK with other brands of resin composite to avoid incomplete cure, entrapment of air bubbles or losing its color adjustment.
- To avoid cross infection, after extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.
- This product is provided non-sterile and is not intended to be sterilized by the user.
- Use OMNICHROMA BLOCKER or OMNICHROMA BLOCKER FLOW under OMNICHROMA FLOW BULK for the direct restoration of extensive Class III or IV cavity. OMNICHROMA FLOW BULK can appear semi-translucent when OMNICHROMA FLOW BULK is not adequately surrounded by a cavity wall. If applied in direct restorations of extensive Class III or IV cavity, the cured OMNICHROMA FLOW BULK may appear dark.
- In the case of repair of porcelain / composite, the application of an opaquer to the fractured surfaces as an initial layer is needed to hide a metal portion if it exists. Then, place OMNICHROMA FLOW BULK onto the fractured surfaces.
- In the case of masking tooth discoloration, apply OMNICHROMA BLOCKER, OMNICHROMA BLOCKER FLOW or a dental color masking material such as an opaquer on the discolored surface to mask the discoloration, followed by the placement of OMNICHROMA FLOW BULK.

## 8. Curing

Light-cure each increment for at least the indicated time (see aforementioned table). Light-cure each layer using Halogen light-curing or LED light-curing unit with a wavelength of 400-500 nm. Hold the light guide tip as close to the restorative as possible during light exposure.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

## 9. Finishing

Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**IMPORTANT NOTE:** The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.

Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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# Packaging waste disposal guide

(According to European Commission Decision 97/129/EC)

## ENGLISH

### Packaging waste disposal guide

|                           |                                                                                          |
|---------------------------|------------------------------------------------------------------------------------------|
| SYRINGE (cap and body)    | Dispose of as industrial waste.                                                          |
| Packaging bag for tip     | LDPE4 (Low density polyethylene):<br>Recycling. Follow local instructions for disposal.  |
| Paperbox and/or partation | PAP21 (Non-corrugated fibreboard):<br>Recycling. Follow local instructions for disposal. |
| IFU and/or guide          | PAP22 (Paper):<br>Recycling. Follow local instructions for disposal.                     |

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